


Respireo

Hospital F

Vented Disposable

- Also available in **vented** version
- Technical features identical to Respireo Hospital Non Vented Disposable
- Fitting instructions cf. recto



Compatibility with ventilators

	Monnal T50	Monnal T60	Monnal T75	Extend	Ventilators (with valve)	Ventilators (with intentional leak)
Respireo Hospital Non Vented	✓	✓	✓	✓	✓	X
Respireo Hospital Vented	X	X	X	X	X	✓

Codes Designation

Size

KM222100	Respireo Hospital F Non Vented Disposable (manufacturer reference 472542*)	S
KM222200	Respireo Hospital F Non Vented Disposable (manufacturer reference 472543*)	M
KM222300	Respireo Hospital F Non Vented Disposable (manufacturer reference 472544*)	L
KM226400	Respireo Hospital F Vented Disposable (manufacturer reference 472539*)	S
KM226500	Respireo Hospital F Vented Disposable (manufacturer reference 472540*)	M
KM226600	Respireo Hospital F Vented Disposable (manufacturer reference 472541*)	L

Respireo
Hospital F
Non Vented Disposable

Respireo
Hospital F
Vented Disposable



Please read user manual
(*) Medical device class IIa CE 0051 manufactured by Air Liquide Medical Systems S.r.l.
Distributed by Air Liquide Medical Systems S.A.

Respireo

Hospital F

Non Vented Disposable



Single-patient use up to 7 days




The solution to benefit from a complete expertise



- Adjustable headgear
- Dual-wall cushion
- Quick-release clips
- Transparent frame
- Reinforced forehead pad
- Adjustable forehead support
- Measurement port

- Provides non-invasive, comfortable ventilation
- Guarantees complete safety for the user
- Saves time by simplifying traceability
- Mask available in 3 sizes
- Sizing-gauge attached on the packaging
- Package includes mask, cushion and headgear

Fitting instruction



- Biocompatible silicone
- Phtalate-free
- Latex-free