

Contact

Air Liquide Medical Systems
Parc de Haute Technologie
6, rue Georges Besse
92182 Antony Cedex, France
Tel.: +33 (0)1 40 96 66 00
Fax : +33 (0)1 40 96 67 00

PS-PRO by Air Liquide Medical Systems

PS-Pro

A self-adapting mode for patient recovery



YN0040411 - V.1.0 - 01.08.2014 - Art direction : Raffaello Fischer

Comfortable

Self-adapting

Continuous

Safe

Comfortable:

Patient is free to breath spontaneously above RRmini

Self-adapting:


From mandatory pressure ventilation to pressure support ventilation

Continuous:

Smooth recovery and maximum tolerance of the treatment

Safe:

3 parameters to secure ventilation: RR, Target Vt, PI max



Mandatory Pressure Cycles


Pressure Support Cycles

Mandatory Pressure Cycles

Pressure Support Cycles



Vt target & Pressure adaptation



PS-Pro progressively decreases pressure level to facilitate safe extubation.

